

# Anti-Glutaminase Rabbit antibody



Catalog Number: 167235

## Product name

Rabbit IgG

Anti-Glutaminase Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, FC

Human

## Dilutions

## Antibody description

WB: 1/2000

Glutaminase Rabbit polyclonal antibody

IHC: 1/50

## Preparation

ICC/IF: 1/50

Antigen: A synthetic peptide of human Glutaminase

FC: 1/50

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type



Western blot detection of Glutaminase-1/GLS1 in Jurkat,SKOV-3,HepG2 cell lysates using Glutaminase-1/GLS1 Rabbit pAb(1:1000 diluted).Predicted band size:73KDa.Observed band size:73,65KDa.