

# DAO antibody



Catalog Number: 107235

## Product name

Monoclonal

DAO antibody

## Ig Type

## Specificity

Mouse IgG1

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

DAO Mouse Monoclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 39 kDa

## Dilutions


Recommended Dilution:

## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of DAO recombinant protein (Accession Number: NM\_001917). Purification method: Protein G purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107235(DAO Antibody) at dilution of 1:1000

## Clonality