

# PRDM5 Polyclonal Antibody



Catalog Number: 166635

## Product name

PRDM5 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to PRDM5

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human PRDM5 (NP\_061169.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations




Western blot - PRDM5 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PRDM5 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - PRDM5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat spleen using PRDM5 antibody at dilution of 1:100



(40x lens).

# PRDM5 Polyclonal Antibody




Catalog Number: 166635

---

Immunohistochemistry - PRDM5 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse liver using PRDM5 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - PRDM5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse spleen using PRDM5 antibody at dilution of



1:100 (40x lens).

Immunofluorescence - PRDM5 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using PRDM5 antibody . Blue: DAPI for nuclear staining.