

# Anti-DR5 (6F6) Mouse antibody



Catalog Number: 168335

## Product name

Anti-DR5 (6F6) Mouse antibody

## Specificity

Human

## Antibody description

DR5 (6F6) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human DR5 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

WB


## Dilutions

WB: 1/1000;Other 1/2000

## Validations



Western blot detection of DR5 in SW480, K562, HCT116 and HeLa cell lysates using DR5 mouse mAb (1:1000 diluted). Predicted band size: 40/48KDa. Observed band size: 40/48KDa.



Western blot detection of DR5 in MCF-7 and T47D cell lysates using DR5 mouse mAb (1:2000 diluted). Predicted band size: 40/48KDa. Observed band size: 40/48KDa.