

REST Polyclonal Antibody



Catalog Number: 163735

Product name

REST Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to REST

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 760-1060 of human REST (NP_001180437.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - REST Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using REST antibody .