

# Anti-FGF2 Rabbit antibody



Catalog Number: 166735

## Product name

Anti-FGF2 Rabbit antibody

## Specificity

Human

## Antibody description

FGF2 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human FGF2

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IP, FC


## Dilutions

WB: 1/500-1/2000

IP: 1/50

FC: 1/50

## Validations



Western blot analysis of FGF2 expression in K562 cell lysate.