

Anti-FGF2 Rabbit antibody



Catalog Number: 166735

Product name

Anti-FGF2 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IP, FC

Antibody description

FGF2 Rabbit polyclonal antibody

Dilutions

WB: 1/500-1/2000

Preparation

Antigen: A synthesized peptide derived from human FGF2


IP: 1/50

FC: 1/50

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Validations



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Western blot analysis of FGF2 expression in K562 cell lysate.