

Anti-LILRB2 / ILT4 / LIR-2 antibody



Catalog Number: 100045

Product name

Anti-LILRB2 / ILT4 / LIR-2 antibody

Immunogen

[Human LILRB2 / ILT4 / LIR-2 \(His Tag\) recombinant protein](#)

Specificity

Human LILRB2 / ILT4 / LIR-2

Antibody description

Mouse monoclonal to LILRB2 / ILT4 / LIR-2

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human LILRB2 / ILT4 / LIR-2 (rh LILRB2 / ILT4 / LIR-2; AAH36827.1; Met1-Val461) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of

water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

FCM

Dilutions

FCM

Validations



LILRB2 / ILT4 / LIR-2 Antibody, Mouse MAb, Flow Cytometry

Flow cytometric analysis of Human LILRB2(CD85d) expression on human whole blood monocytes. Cells were stained with FITC-conjugated anti-Human LILRB2(CD85d). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable monocytes.