

# DOCK6 antibody



Catalog Number: 110045

## Product name

DOCK6 antibody

## Immunogen

[Human DOCK6 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

DOCK6 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 230 kDa

## Preparation

This antibody was obtained by immunization of DOCK6 recombinant protein (Accession Number: NM\_020812). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:110045(Dock6 Antibody) at dilution of 1:300