

Anti-JAG1/Jagged 1(CD339) antibody



Catalog Number: 101045

Product name

Anti-JAG1/Jagged 1(CD339) antibody

Immunogen

[Human JAG1/Jagged 1\(CD339\) \(His Tag\) recombinant protein](#)

Specificity

Human JAG1 / Jagged 1

Antibody description

Mouse monoclonal to JAG1/Jagged 1(CD339)

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human JAG1 / Jagged 1 (rh JAG1 / Jagged 1; NP_000205.1; Met1-Ser1046) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !

Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG2a


Applications

FCM

Dilutions

FCM

Validations



JAG1 / Jagged 1 Antibody (FITC), Mouse MAb,
Flow cytometric

Profile of anti-JAG1(CD339) reactivity on HepG2 cells analyzed by flow cytometry.