

H7N9 HA (His Tag) recombinant protein



Catalog Number: 501045

General Information

Protein Construction

A DNA sequence encoding the influenza A virus (A/Zhejiang/DTID-ZJU10/2013(H7N9)) hemagglutinin (AHA11500.1) (Met1-Val524), termed as HA, was fused with a polyhistidine tag at the C-terminus.

Organism

H7N9

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 19

Molecular Mass

The secreted recombinant HA subunit of influenza A virus (A/Zhejiang/DTID-ZJU10/2013(H7N9))

comprises 517 amino acids and has a predicted molecular mass of 57.6 kDa. It migrates as an approximately 63.8 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H7N9 (A/Zhejiang/ DTID-ZJU10/2013) Hemagglutinin / HA Protein (His Tag) SDS-PAGE