

KIAA0182 antibody



Catalog Number: 112045

Product name

KIAA0182 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

KIAA0182 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IHC detected in human testis tissue, human kidney tissue. Observed molecular weight by Western-blot: 136 kDa

Preparation

This antibody was obtained by immunization of KIAA0182 recombinant protein (Accession Number: XM_047433833). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112045(KIAA0182 Antibody) at dilution of



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:112045(KIAA0182 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:112045(KIAA0182 Antibody) at dilution of 1:50

KIAA0182 antibody



Catalog Number: 112045

(under 40x lens)