

# KIAA0182 antibody



Catalog Number: 112045

## Product name

KIAA0182 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

KIAA0182 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IHC detected in human testis tissue, human kidney tissue. Observed molecular weight by Western-blot: 136 kDa

## Preparation

This antibody was obtained by immunization of KIAA0182 recombinant protein (Accession Number: XM\_047433833). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112045(KIAA0182 Antibody) at dilution of



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:112045(KIAA0182 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:112045(KIAA0182 Antibody) at dilution of 1:50

# KIAA0182 antibody



Catalog Number: 112045

---

(under 40x lens)