

# SAG antibody

Catalog Number: 115045

---

## Product name

SAG antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SAG Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of SAG recombinant protein (Accession Number: NM\_000541). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:115045(SAG Antibody) at dilution of 1:600