

Factor XII antibody



Catalog Number: 110445

Product name

Factor XII antibody

Recommended Dilution:

WB: 1:200-1:2000

Immunogen

IHC: 1:20-1:200

[Human F12 Recombinant protein \(His tag\)](#)


Validations

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Factor XII Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human colon cancer tissue. Positive WB detected in Jurkat cells, HepG2 cells, mouse liver tissue. Observed molecular weight by Western-blot: 28 kDa




Preparation

This antibody was obtained by immunization of Factor XII recombinant protein (Accession Number: BC012390). Purification method: Antigen affinity purified.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110445(F12 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Immunohistochemical of paraffin-embedded human liver using Catalog No:110445(F12 antibody) at dilution of 1:50 (under 10x lens)



Dilutions

Immunohistochemical of paraffin-embedded

Factor XII antibody



Catalog Number: 110445

human liver using Catalog No:110445(F12

antibody) at dilution of 1:50 (under 40x lens)