

ERRFI1 Polyclonal Antibody



Catalog Number: 161445

Product name

ERRFI1 Polyclonal Antibody

Specificity

Rat

Antibody description

Polyclonal antibody to ERRFI1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 223-462 of human ERRFI1 (NP_061821.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - ERRFI1 Polyclonal Antibody

Western blot analysis of extracts of rat lung, using ERRFI1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 5min.