

# Anti-PPP1R15B antibody



Catalog Number: 176445

## Product name

Anti-PPP1R15B antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to PPP1R15B

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human ppp1r15b. aa range:411-460

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

## Validations





Fig1:: Western Blot analysis of NIH-3T3 cells using PPP1R15B Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at



1:20000

Fig2:: Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100