

p504S,AMACR antibody



Catalog Number: 107445

Product name

p504S,AMACR antibody

Specificity

human, mouse; other species not tested.

Antibody description

p504S,AMACR Mouse Monoclonal antibody.
Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 42kd

Preparation

This antibody was obtained by immunization of p504S,AMACR recombinant protein (Accession Number: NM_203382). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107445(p504S,AMACR antibody) at dilution of 1:1000