

# GALNS antibody

Catalog Number: 110845

## Product name

GALNS antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GALNS Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue. Positive WB detected in human liver tissue. Observed molecular weight by Western-blot: 50-55kd

## Preparation

This antibody was obtained by immunization of GALNS recombinant protein (Accession Number: NM\_000512). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110845(GALNS antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human kidney using Catalog No:110845(GALNS antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:110845(GALNS antibody) at dilution of 1:50 (under 40x lens)