

IREB2 Polyclonal Antibody



Catalog Number: 164845

Product name

IREB2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to IREB2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human IREB2 (NP_004127.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:20 - 1:50

Validations



Immunofluorescence - IREB2 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using IREB2 antibody . Blue: DAPI for nuclear staining.