

Anti-Rab5 Rabbit antibody



Catalog Number: 166845

Product name

Anti-Rab5 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Rab5 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Rab5

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/50

Validations



Western blot detection of Rab5 in K562, C6, 3T3, HeLa cell lysates using Rab5 Rabbit pAb(1:1000 diluted). Predicted band size: 24kDa. Observed band size: 24kDa.