

# Anti-Cytokeratin 16 Rabbit antibody



Catalog Number: 168845

## Product name

Anti-Cytokeratin 16 Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IHC-P

## Antibody description

Cytokeratin 16 Rabbit monoclonal antibody

## Dilutions

1:2000

## Preparation


Antigen: Recombinant protein of human Cytokeratin 16

1:50

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Cytokeratin 16 in HeLa cell lysates using Cytokeratin 16 Rabbit mAb(1:1000 diluted).Predicted band size:51kDa.Observed band size:51kDa.

## Clonality

Monoclonal