

DECR2 antibody



Catalog Number: 109845

Product name

DECR2 antibody

Immunogen

[Human DECR2 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

DECR2 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of DECR2 recombinant protein (Accession Number: BC011968). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:109845(DECR2 Antibody) at dilution of 1:1000