

DYRK4 antibody



Catalog Number: 110145

Product name

DYRK4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DYRK4 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, human brain tissue, Jurkat cells, mouse testis tissue, U-937 cells. Positive IHC detected in human lymphoma tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of DYRK4 recombinant protein (Accession Number: NM_003845). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110145(DYRK4 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:110145(DYRK4 antibody) at dilution of 1:100 (under 10x lens)