

NFE2L1 antibody



Catalog Number: 113145

Product name

NFE2L1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

NFE2L1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of NFE2L1 recombinant protein (Accession Number: NM_001330262). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113145(NFE2L1 Antibody) at dilution of 1:600