

Anti-CALML5 antibody



Catalog Number: 106145

Product name

Anti-CALML5 antibody

Immunogen

[Human CALML5 \(His & GST Tag\) recombinant protein](#)

Specificity

Human CALML5

Antibody description

Rabbit monoclonal to CALML5

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human CALML5 (rh CALML5; AAH39172.1; Met1-Glu146).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 0.5-5 µg/mL

Validations



CALML5 Antibody, Rabbit MAb,
Immunohistochemistry

Immunochemical staining of human CALML5 in human skin with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to epithelium.