

Cosibelimab reference antibody



Catalog Number: 177145

Product name

Cosibelimab reference antibody

Specificity

Human B7-H1 / PD-L1 / CD274

Antibody description

Anti-B7-H1 / PD-L1 / CD274 Antibody
(Cosibelimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations



SDS-PAGE

Anti-B7-H1 / PD-L1 / CD274 Antibody (cosibelimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-B7-H1 / PD-L1 / CD274 Antibody (cosibelimab) is 99.9%, determined by SEC-HPLC.