

FLII Polyclonal Antibody



Catalog Number: 162545

Product name

FLII Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to FLII

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 400-660 of human FLII (NP_002009.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:200 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - FLII Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using FLII antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5S.



Immunohistochemistry - FLII Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using FLII antibody at dilution of 1:100 (40x




FLII Polyclonal Antibody



Catalog Number: 162545

Immunohistochemistry - FLII Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using FLII antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - FLII Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using FLII antibody at dilution of 1:100 (40x lens).