

RAP1B-Specific antibody



Catalog Number: 114545

Product name

RAP1B-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAP1B-Specific Rabbit Polyclonal antibody.
Positive WB detected in human brain tissue, A549, mouse brain tissue. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of RAP1B-Specific recombinant protein (Accession Number: NM_001010942). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114545(RAP1B antibody) at dilution of 1:300



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:114545(RAP1B antibody) at dilution of 1:300. (Data provided by Angran Biotech (www.miRNAlab.com)).