

Phospho-PEA15-S104 pAb



Catalog Number: 166545

Product name

Phospho-PEA15-S104 pAb

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Phospho-PEA15-S104

Preparation

Antigen: A synthetic phosphorylated peptide around S104 of human PEA15 (NP_003759.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations




Immunohistochemistry - Phospho-PEA15-S104 pAb

Immunohistochemistry of paraffin-embedded rat lung using Phospho-PEA15-S104 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - Phospho-PEA15-S104 pAb

Immunohistochemistry of paraffin-embedded rat testis using Phospho-PEA15-S104 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - Phospho-PEA15-S104

Phospho-PEA15-S104 pAb

Catalog Number: 166545


pAb

Immunohistochemistry of paraffin-embedded rat brain using Phospho-PEA15-S104 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - Phospho-PEA15-S104 pAb

Immunohistochemistry of paraffin-embedded human tonsil using Phospho-PEA15-S104 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - Phospho-PEA15-S104 pAb

Immunohistochemistry of paraffin-embedded human lung cancer using Phospho-PEA15-S104 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - Phospho-PEA15-S104 pAb

Immunohistochemistry of paraffin-embedded human appendix using Phospho-PEA15-S104 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - Phospho-PEA15-S104 pAb

Immunohistochemistry of paraffin-embedded human gastric cancer using Phospho-PEA15-S104 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - Phospho-PEA15-S104


Phospho-PEA15-S104 pAb



Catalog Number: 166545

pAb


Immunohistochemistry of paraffin-embedded mouse lung using Phospho-PEA15-S104 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - Phospho-PEA15-S104 pAb

Immunohistochemistry of paraffin-embedded

mouse testis using Phospho-PEA15-S104 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - Phospho-PEA15-S104 pAb

Immunohistochemistry of paraffin-embedded mouse brain using Phospho-PEA15-S104 antibody at dilution of 1:100 (40x lens).