

Anti-E Cadherin Rabbit antibody



Catalog Number: 167545

Product name

Anti-E Cadherin Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Mouse, Rat

Applications

WB, IP

Antibody description

E Cadherin Rabbit polyclonal antibody

Dilutions

WB: 1/1000

IP: 1/20


Preparation

Antigen: Recombinant protein of mouse E Cadherin

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of E Cadherin in Mouse brain lysates using E Cadherin Rabbit pAb(1:1000 diluted).Predicted band size:98(cleavages),135KDa.Observed band size: 80-120(cleavages),135KDa.

Clonality

Polyclonal