

Anti-VEGFD Rabbit antibody



Catalog Number: 168945

Product name

Anti-VEGFD Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

VEGFD Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human VEGFD

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


1:1000-1:5000

1:50


1:50

1:20-1:100

Validations



Immunocytochemistry of VEGFD(green) in HeLa cells using VEGFD Rabbit mAb at dilution 1/50,



Western blot detection of VEGFD in Rat Brain,C6,HeLa cell lysates using VEGFD Rabbit mAb(1:1000 diluted).Predicted band size:40kDa.Observed band size:40kDa.