

DIRAS2 antibody



Catalog Number: 109945

Product name

DIRAS2 antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

[Human DIRAS2 Recombinant protein \(GST tag\)](#)

IP: 1:200-1:2000

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

DIRAS2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells. Positive IP detected in mouse brain tissue. Positive FC detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 22 kDa




Preparation

This antibody was obtained by immunization of DIRAS2 recombinant protein (Accession Number: NM_017594). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109945(DIRAS2 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunofluorescent analysis of HeLa cells, using DIRAS2 antibody Catalog No:109945 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

Ig Type

Rabbit IgG



Applications

ELISA, WB, IF, FC, IP

Dilutions

DIRAS2 antibody

Catalog Number: 109945



IP Result of anti-DIRAS2 (IP:Catalog No:109945, 4ug; Detection:Catalog No:109945 1:500) with mouse brain tissue lysate 4000ug.

1×10^6 HeLa cells were stained with 0.05ug DIRAS2 antibody (Catalog No:109945, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488 - Goat anti-Rabbit IgG with dilution 1:100.