

Anti-SNX6 antibody



Catalog Number: 142245

Product name

Anti-SNX6 antibody

Specificity

Human SNX6

Antibody description

Rabbit Polyclonal to Human SNX6

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human SNX6, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 10-20 µg/ml

Validations



Items	Lanes	A
Sample (whole cell lysate)		HepG2
Sample Volume (µg/lane)		60
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		10-20 µg/ml
Secondary Antibody		Dylight 800 labeled Antibody to Rabbit IgG [H+], at 1:5000 dilution.
Explanation		Developed using Odyssey imaging system. Predicted band size : 47 kDa Observed band size : 47 kDa Additional bands at : 34 kDa (We are unsure as to the identity of these extra bands.)

Human SNX6 Western blot (WB) 14561