

# MT-CO3 Polyclonal Antibody



Catalog Number: 166245

## Product name

MT-CO3 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to MT-CO3

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse MT-CO3 (NP\_904334.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



## Western blot - MT-CO3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MT-CO3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.