

Anti-PRKACG Rabbit antibody



Catalog Number: 167245

Product name

Anti-PRKACG Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PRKACG Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human PRKACG

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, FC, IP


Dilutions

WB: 1/5000-1/10000

FC: 1/20-1/100

IP: 1/20-1/50

Validations



Western blot detection of PRKACG in MCF7,C6,HeLa cell lysates using PRKACG Rabbit pAb(1:1000 diluted).Predicted band size:40KDa.Observed band size:40KDa.