

Anti-PRTFDC1 antibody



Catalog Number: 100645

Product name

Anti-PRTFDC1 antibody

Immunogen

[Human PRTFDC1 \(His Tag\) recombinant protein](#)

Specificity

Human PRTFDC1

Antibody description

Mouse monoclonal to PRTFDC1

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human PRTFDC1 (rh PRTFDC1; Q9NRG1-1; Met1-Val225). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% treha

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG


Applications

WB

Dilutions

WB: 20-30 µg/ml

Validations



Lanes	A
Items	
Sample (Membrane Lysate)	293T
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	20-30µg/ml
Secondary Antibody	Dylight 800 labeled Antibody to Mouse IgG (H+L), at 1:2500 dilution.
	Developed using Odyssey imaging system.
Explanation	Predicted band size : 26 kDa Observed band size : 26 kDa

PRTFDC1 Antibody, Mouse MAb, Western blot