

# Human PFN4 / Profilin 4 (His Tag) recombinant protein



Catalog Number: 500645

## General Information

### Gene Name Synonym

Profilin IV

### Protein Construction

A DNA sequence encoding the human PFN4 (NP\_955378.1) (Met1-Ser129) was expressed with a polyhistidine tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

His

### Molecular Mass

The recombinant human PFN4 consists of 144

amino acids and predicts a molecular mass of 16.2 KDa. It migrates as an approximately 14 KDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, 200mM Arg, pH 7.4.  
1. 5 % trehalose and mannitol are added as protectants before lyophilization.  
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human PFN4 / Profilin 4 Protein (His Tag) SDS-PAGE