

PDZK1 antibody



Catalog Number: 113645

Product name

PDZK1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PDZK1 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in mouse kidney tissue, HEK-293 cells. Observed molecular weight by Western-blot: 63-70 kDa

Preparation

This antibody was obtained by immunization of PDZK1 recombinant protein (Accession Number: XM_017001464). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113645(PDZK1 antibody) at dilution of 1:500



IP Result of anti-PDZK1 (IP:Catalog No:113645, 3ug; Detection:Catalog No:113645 1:500) with HEK-293 cells lysate 3200ug.