

Product name

ARSB Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to ARSB

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 404-533 of human ARSB (NP_000037.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Immunohistochemistry - ARSB Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using ARSB antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ARSB Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using ARSB antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - ARSB Polyclonal Antibody

ARSB Polyclonal Antibody



Catalog Number: 163345

Immunohistochemistry of paraffin-embedded mouse liver using ARSB antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ARSB Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using ARSB antibody at dilution of 1:100 (40x lens).