

Anti-GTF3C2 antibody



Catalog Number: 176345

Product name

Polyclonal

Anti-GTF3C2 antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IHC-P

Antibody description

Rabbit polyclonal antibody to GTF3C2

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tf3c2. aa range:71-120

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1:: Western Blot analysis of various cells using TFIIC110 Polyclonal Antibody