

TRAPPC4 antibody



Catalog Number: 116345

Product name

TRAPPC4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TRAPPC4 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 24kd

Preparation

This antibody was obtained by immunization of TRAPPC4 recombinant protein (Accession Number: NM_016146). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116345(TRAPPC4 antibody) at dilution of 1:400