

Human UROS (His tag) recombinant protein



Catalog Number: 519345

General Information

Gene Name Synonym

Uroporphyrinogen-III Synthase, UROIIIIS, UROS, Hydroxymethylbilane Hydrolyase [Cyclizing], Uroporphyrinogen-III Cosynthase, UROS

Protein Construction

A DNA sequence encoding the Human UROS Met1-Cys265 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 10% Glycerol, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human UROS (His tag) recombinant protein