

# IL-4 antibody

Catalog Number: 111745

## Product name

Polyclonal

IL-4 antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

IL-4 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of IL-4 recombinant protein (Accession Number: NM\_000589). Purification method: Antigen affinity purified.

## Validations

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111745(IL4 Antibody) at dilution of 1:700

## Clonality