

Anti-Bcl-2 Rabbit antibody



Catalog Number: 168745

Product name

Anti-Bcl-2 Rabbit antibody

Specificity

Human

Antibody description

Bcl-2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Bcl-2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

1:1000


1:20-1:500

1:100

1:50

1:30

Validations



Western blot detection of Bcl-2 in Jurkat, HeLa cell lysates using Bcl-2 Rabbit mAb(1:1000 diluted). Predicted band size: 26kDa. Observed band size: 26kDa.