

DOK7 antibody

Catalog Number: 110055

Product name

DOK7 antibody

Specificity

Human; other species not tested.

Antibody description

DOK7 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, MCF-7 cells. Positive IHC detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 64 kDa

Preparation

This antibody was obtained by immunization of DOK7 recombinant protein (Accession Number: NM_173660). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110055(DOK7 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:110055(DOK7 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:110055(DOK7 Antibody) at dilution of 1:50 (under 10x lens)

DOK7 antibody



Catalog Number: 110055

No:110055(DOK7 Antibody) at dilution of 1:50 (under 40x lens)