

ZIC4-Specific antibody



Catalog Number: 117055

Product name

ZIC4-Specific antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZIC4-Specific Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 37kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_032153). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117055(ZIC4-Specific antibody) at dilution of 1:500