

Anti-TNS2 antibody



Catalog Number: 176455

Product name

Anti-TNS2 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to TNS2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tensin-2. aa range:439-488

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western blot analysis of A549 3T3 using Tensin-2 antibody. Antibody was diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000