

# TXNL1 antibody



Catalog Number: 116455

## Product name

TXNL1 antibody

## Immunogen

[Human TXNL1 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TXNL1 Rabbit Polyclonal antibody. Positive IP detected in HL-60 cells. Positive WB detected in K-562 cells, HL-60 cells. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of TXNL1 recombinant protein (Accession Number: XM\_017026093). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116455(TXNL1 antibody) at dilution of 1:1000



IP Result of anti-TXNL1 (IP:Catalog No:116455, 4ug; Detection:Catalog No:116455 1:300) with HL-60 cells lysate 960ug.