

Anti-Cardiac Troponin I (2D10) Mouse antibody



Catalog Number: 168455

Product name

Mouse IgG2a

Anti-Cardiac Troponin I (2D10) Mouse antibody

Applications

Specificity

WB, ELISA

Human, Mouse, Rat

Dilutions

Antibody description

ELISA: 1/10000

Cardiac Troponin I (2D10) Mouse monoclonal antibody

;WB: 1/1000

Preparation

Validations

Antigen: Synthetic peptides corresponding to the sequence of human Cardiac Troponin I.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.


Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot analysis of extracts from Rat and mouse heart tissue using Troponin I (168455, dilution 1:1000) Mouse mAb. Predicted band size: 28kDa. Observed band size: 28kDa.